
CURING LIGHT

LED-2LIGHT

fr LAMPE LED À PHOTOPOLYMÉRISER SANS FIL
en .. WIRELESS LED LIGHT CURING UNIT
it..... LAMPADA LED PER POLIMERIZZAZIONE
SENZA FILO
es... LÁMPARA LED SIN CABLE
pt ... FOTOPOLIMERIZADOR LED SEM FIOS
de .. DRAHTLOSE LED POLYMERISATIONSLAMPE
nl.... DRAADLOZE LED LICHTUITHARDINGSLAMP

fr : Guide d'utilisation.....	03 - 10
en: User's guide.....	11 - 18
it : Manuale d'uso	19 - 26
es : Manual de usuario	27 - 34
pt : Guia do utilizador.....	35 - 42
de: Benutzerhandbuch.....	43 - 50
nl : Bedienungsanleitung.....	51 - 58

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Description :

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PRODUCT OVERVIEW

Description:

The curing light uses the principle of radiation to cure light-curing resin by briefly irradiating it; it is intended for use in hospitals and dental practices.

This product comes with a protective screen, an light guide, a curing handpiece, a charging station, a mains adapter and 100 single-use protective sleeves.

List of accessories: Protective screen, single-use protective sleeves, battery, light guide, battery cover, mains adapter, charging station.

Target population:

Adults and children

Intended use:

This curing light is intended for use by qualified dental professionals for the light-curing of dental resins and composite materials.


CONTRAINDICATIONS

Patients with heart conditions, pregnant women, children and people allergic to blue-violet light should exercise caution when using this device.

TECHNICAL SPECIFICATIONS

Light guide	10 mm (diameter)	
Wavelength	385-515 nm	
Modes	Soft	700 to 900 mW/cm ² for 5, 10, 15 or 20 seconds
	Ramp (Progressive)	1,000 to 1,200 mW/cm ² for 5, 10, 15 or 20 seconds
	Normal	1,000 to 1,200 mW/cm ² for 5, 10, 15 or 20 seconds
	Ortho	2,300 to 2,500 mW/cm ² for 3 seconds at 10 pulses per cycle, or 5 seconds at 10 pulses per cycle with a 1-second interval
	High (High power)	2,300 to 2,500 mW/cm ² for 1, 3 or 5 seconds
	Turbo	2,700 to 3,000 mW/cm ² for 1, 2 or 3 seconds
Electrical specifications	Mains adapter input: 100–240 V AC, 50 Hz/60 Hz, 0.4 A max. Mains adapter output: 5 V DC, 1.0 A Internal battery: 3.7 V, 1,900 mAh	
Software version	V1.0.0	
Environmental factors	Operation: 5 °C to 40 °C / 30% to 75% relative humidity / 70 kPa to 106 kPa atmospheric pressure Transport and storage: -20 °C to +55 °C / 10% to 93% relative humidity / 70 kPa to 106 kPa atmospheric pressure	
Equipment safety classification	1. Operating mode: Non-continuous operation: the equipment can only operate continuously for a maximum of 20 second 2. Type of protection against electric shock: Class II 3. Degree of protection against electric shock: Type B applied part 4. Applied part of the equipment: light guide 5. Degree of protection against water ingress: IPX0 6. Safety level for use in explosive atmospheres or in the presence of flammable anaesthetic gases (air, oxygen, nitrous oxide): Use prohibited	

INSTALLATION

1. Charging:

Plug the power cable directly into the curing light's charging station. Place the main unit on the charging base; the curing light will then begin to charge. The charge level is displayed on the screen and scrolls until charging is complete. Low battery warning: The battery symbol on the display will start to flash; we recommend that you recharge it in good time.

2. Curing light handpiece:

- Insert the light guide into the curing light handpiece BY TWISTING IT. Ensure that the light guide is securely in place in the curing light handpiece.
- The battery is not fully charged prior to dispatch. Please charge the battery for at least 4 hours before first use.
- The curing light's circuitry automatically compensates for voltage fluctuations during use to ensure a constant light output.
- Leaving the curing light on charge at all times will not affect the battery's lifespan, as it is equipped with a built-in safety feature that stops charging as soon as it is fully charged.
- Light intensity must be checked regularly using a lux meter.

OPERATION

- Switch on the curing light by pressing the ON/OFF switch ONCE.
- The curing light is factory-preset to 'Normal' mode for 10 seconds.
- Press the MODE switch to change modes. Press the DURATION switch several times to select the desired duration for the chosen mode.
- The start and end of each programme, as well as any change in duration or mode, will be signalled by a 1 "beep".
- The curing light switches off automatically if it remains inactive for 2 minutes. It retains the last selected programme in its memory once it is restarted.
- Slip the single-use protective sleeve over the curing light so that it is completely covered to prevent cross-contamination.
- Place the protective screen over the light guide of the curing light until it is securely attached. The curing light is now ready for use.
- Press the ON/OFF switch to switch on the light-curing unit if none of the indicator lights are on. Once preparations for light curing are complete, press the ON/OFF switch to start using the light-curing unit.
- Once the light-curing process is complete, remove and dispose of the single-use protective sleeve after each patient.

TROUBLESHOOTING

Fault	Possible cause	Solutions
No response	1. Battery flat 2. Faulty battery 3. Battery protection	1. Recharge the battery or replace it with a new one 2. Replace the battery with a new one 3. Activate it by placing it on the base
Insufficient light intensity	Resin remains on the surface of the light guide	Remove the resin
Does not charge when connected to the adapter	1. The adapter is not properly plugged in 2. Faulty or incompatible adapter 3. Dirt on the charging base connector	1. Plug it in again 2. Replace the adapter 3. Clean the internal connector of the charging base with alcohol
Reduced battery life after a full charge	Decreased battery capacity	Replace the battery with a new one
The mode indicator is flashing when the battery icon flashes on the main unit's screen	Low voltage	Returns to normal after 15 minutes of charging

If you have carried out these steps but the device is still not working normally, please contact your local dealer or the manufacturer.

This information is provided for guidance only; please always follow the instructions provided by the manufacturers of your resin and composite material.

- Numerous factors (distance and angle between the tip of the light guide and the composite material, surface, shade, thickness and type, etc.) affect the curing times of the composite material. It is recommended to apply and light-cure only 2 mm of composite material at a time.
- Using a protective sleeve reduces the light output intensity by 5–10 per cent.
- The light guide of the lamp should be positioned approximately 1–2 mm from the resin or composite material and oriented perpendicular to the tooth surface.
- High-intensity curing lights generate more heat, which is a direct indication of the power of the light emitted. When using the curing light for more than 10 seconds, do not hold the tip of the light guide in the same position. It is always recommended to use a dental syringe to direct a stream of cold air onto the affected area during prolonged periods of light-curing. This will help to keep the temperature of the composite material and the surrounding area lower, particularly when using a rubber dam, in order to prevent the dam from rupturing due to heat build-up.
- Do not expose soft tissues to the curing light for more than 2 seconds, as this may cause burns.

CLEANING AND DISINFECTION

Switch off the curing light and unplug the power cord from the adapter. The curing light handpiece, light guide and protective screen may be cleaned with clean water or a neutral sterilising solution, for cold disinfection purposes only. Do not leave them to soak.

⚠ Warnings

- For use exclusively by trained dental professionals, in accordance with its intended use.
- This device must be used in strict accordance with these INSTRUCTIONS. The manufacturer accepts no liability in the event of failure to comply with the INSTRUCTIONS or use of the device for any other purpose.
- The light emitted by this type of device may be hazardous and MUST NOT UNDER ANY CIRCUMSTANCES be directed towards the eyes. The light emitted by this device must be directed exclusively towards the area to be treated within the oral cavity. Always use protective screens and safety goggles whilst the curing light is in operation.
- Do not use the device near a heat source. Do not use solvents, detergents or flammable products to clean the device or immerse it in such substances, as this could damage it or cause a short circuit.
- Ensure that no liquid enters the curing light or the adaptor.
- Never modify the device or any of its components. Any modification could compromise its safety and effectiveness.
- The curing lamp cannot be repaired in the field. Do not dismantle this product. Unauthorised repair or tampering with internal parts may result in serious injury and VOIDS the manufacturer's warranty.
- Check the device before each use. Do not use the device if it shows any signs of damage. Prolonged use of a damaged device may cause injury or produce incorrect results.
- Do not place the device in a position where it would be difficult to disconnect it from the mains supply.
- Charge the device for at least 2 hours before use if it has not been used for more than a month.
- The curing light is fitted with an overheating protection system. It can operate continuously for 200 seconds; for example, if you use the curing light 10 times in succession in 20-second bursts, it will then switch to overheating protection mode. Only after 2 minutes of inactivity will it be able to resume continuous operation for 200 seconds.
- Only the original charger and adaptor may be used, as chargers and adaptors from other brands may damage the circuitry. Only the lithium battery specified or supplied by the manufacturer may be used to ensure the device functions correctly.

TRANSPORT, STORAGE AND MAINTENANCE

Transport

During transport, avoid knocks and vibrations; handle the device with care and do not turn it upside down. Do not transport it together with dangerous goods. Avoid exposure to sunlight, rain or snow during transport.

Storage

Do not store it with toxic, corrosive, flammable or explosive materials. The product must be stored in an environment where the relative humidity does not exceed 10% to 93%, the atmospheric pressure is between 70 kPa and 106 kPa, and the temperature is between -20 °C and +55 °C.

Maintenance

This equipment must be handled with care, kept away from any source of seismic vibrations, and must be installed or stored in a cool, dry and well-ventilated place. When the device is not in use, switch it off using the power switch and unplug the mains plug; if it is not going to be used for a prolonged period, recharge the main unit every 3 months. To clean or disinfect it, wipe the surface with clean water or a disinfectant; do not immerse it in any liquid. After use, check that no resin remains on the surface of the lamp, so as not to compromise the service life of the light guide or the effectiveness of the polymerisation.

EMC DECLARATION OF CONFORMITY

Instructions for Use

The ME EQUIPMENT or ME SYSTEM is intended for use in hospitals or dental practices.

Warning: Do not place this equipment near operating HF surgical devices or within the RF-shielded room of an ME system intended for magnetic resonance imaging, where the intensity of electromagnetic interference is high.

Warning: This equipment should not be used in the immediate vicinity of, or stacked on top of, other equipment, as this may result in malfunction. If such use is necessary, this equipment and the other equipment should be monitored to ensure they are operating normally.

Warning: The use of accessories, transducers and cables other than those specified or supplied by the manufacturer of this equipment may result in increased electromagnetic emissions or reduced electromagnetic immunity of this equipment, and may lead to malfunction.

Warning: Portable RF communication equipment (including peripherals such as antenna cables and external antennas) must not be used within 30 cm (12 inches) of any part of the equipment, including cables specified by the manufacturer. Failure to do so may result in a deterioration in the performance of this equipment.

Note: The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). When used in a residential environment (where Class B of the CISPR 11 standard is normally required), this equipment may not provide sufficient protection for radio frequency communication services. The user may need to take mitigation measures, such as moving or reorienting the equipment.

List of all cables

N°	Name	Length	Shielded or unshielded	Detachable or non-detachable
1	Mains adapter output cable	1.2 m	No	No

Replaceable accessories

N°	Name	Model / Specifications	Connection type
1	Mains adapter	5 V/1 A	Connector
2	Battery	ICR 18500	Plug-in connector

Essential performance characteristics

The curing light has no vital functions nor any diagnostic functions related to the maintenance of vital functions.

The following functions are observed:

- Curing of dental resins and composite materials
- Continuous light emission from the handpiece

If the product is faulty or impaired due to electromagnetic interference, the user must stop using it immediately to avoid any errors caused by a malfunction or a deterioration in its performance. In this case, the user must eliminate the source of interference or adjust the orientation or position of the product so that it can operate normally.

Technical description

Portable and mobile RF communication equipment may affect the performance of the device; strong electromagnetic interference should be avoided during use, and the device should not be placed near mobile phones, microwave ovens, etc. The equipment should not be used in the immediate vicinity of, or stacked on top of, other equipment, as this may result in malfunction. If such use is necessary, this equipment and the other equipment should be monitored to ensure they are operating normally.

With the exception of cables sold by manufacturers as spare parts for internal components, the use of accessories and cables other than those specified or supplied by the manufacturer may result in increased electromagnetic emissions or reduced electromagnetic immunity of this equipment, and may cause malfunction.

The use of accessories, transducers and cables other than those specified or supplied by the manufacturer with this equipment could result in increased electromagnetic emissions or reduced electromagnetic immunity of this equipment, and cause it to malfunction.

Table 1

Manufacturer's recommendations and declaration – Electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class A
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable

Table 2

Manufacturer's recommendations and declaration – Electromagnetic immunity		
Immunity test	Test level IEC 60601-1-2	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air
Fast electrical transients in bursts – IEC 61000-4-4	±2 kV on power supply lines ±1 kV on signal inputs/outputs, 100 kHz repetition frequency	±2 kV on power supply lines Not applicable 100 kHz repetition frequency
Surge voltages – IEC 61000-4-5	±0.5 kV, ±1 kV in differential mode ±0.5 kV, ±1 kV, ±2 kV in common-mode	±0.5 kV, ±1 kV in differential mode; Not applicable
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	0 % UT ; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT ; 1 cycle at 70% UT; 25/30 cycles; Single-phase: at 0°. 0 % UT ; 250/300 cycles	0 % UT ; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT ; 1 cycle at 70% UT; 25/30 cycles; Single-phase: at 0°. 0 % UT ; 250/300 cycles
Power-frequency magnetic field IEC 61000-4-8	30 A/m 50 Hz/60 Hz	30 A/m 50 Hz/60 Hz
Conducted RF IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in the ISM bands between 0.15 MHz and 80 MHz, 80 % amplitude modulation (AM) at 1 kHz	3 V 0.15 MHz – 80 MHz 6 V in the ISM bands between 0.15 MHz and 80 MHz, 80 per cent amplitude modulation (AM) at 1 kHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80 % amplitude modulation (AM) at 1 kHz	3 V/m 80 MHz – 2.7 GHz 80 % amplitude modulation (AM) at 1 kHz
NOTE: UT corresponds to the mains AC voltage prior to the application of the test level.		

Table 3

Recommendations and declaration – Electromagnetic immunity						
Radiated RF IEC 61000-4-3 (Test specifications for immunity of the enclosure port regarding RF wireless communication equipment)	Test frequency (MHz)	Band (MHz)	Service	Modulation	IEC 60601-1-2 standard Test level (V/m)	Compliance level (V/m)
	385	380–390	TETRA 400	18 Hz pulse modulation	27	27
	450	430 - 470	GMRS 460, FRS 460	FM, ±5 kHz deviation, 1 kHz sine wave	28	28
	710	704 - 787	LTE Bands 13 and 17	217 Hz pulse modulation	9	9
	745					
	780					
	810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	18 Hz pulse modulation	28	28
	870					
	930					
	1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE bands 1, 3, 4, 25; UMTS	217 Hz pulse modulation	28	28
	1845					
	1970					
	2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE band 7	217 Hz pulse modulation	28	28
	5240	5100 - 5800	WLAN 802.11 a/n	217 Hz pulse modulation	9	9
	5500					
	5785					

Table 4

Guidelines and manufacturer's declaration – Electromagnetic immunity				
Radiated RF IEC 61000-4-39 (Test specifications for immunity of the enclosure port regarding proximity magnetic fields)	Test frequency	Modulation	Test level in accordance with IEC 60601-1-2 (A/m)	Compliance level (A/m)
	30 kHz	Continuous wave	8	8
	134.2 kHz	Pulse 2.1 kHz modulation	65	65
	13.56 kHz	Pulse 50 kHz modulation	7.5	7.5

ENVIRONMENTAL PROTECTION

At the end of the device's service life, please dispose of it in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive and your country's regulations on the disposal of medical waste.

DISCLAIMER

The manufacturer reserves the right to modify the product at any time without prior notice. The illustrations are provided for guidance only. The right of final interpretation rests with **Guilin Mecco Medical Instrument Co., Ltd.** The industrial design, internal structure, etc., are the subject of several patent applications filed by **Mecco**; any copying or counterfeiting will result in legal action.

WARRANTY

The curing light is guaranteed to the original purchaser against any defects in workmanship and materials, subject to correct installation, use and maintenance, for a period of 12 months from the date of purchase

AFTER-SALES SERVICE FORM

Customer name	
Address	
Postcode	Tel.
Email	
Date of purchase	Date of manufacture
Distributor	
Model	Product reference

EXPLANATION OF SYMBOLS

Warning

Date of manufacture

Class II equipment

Screw in/out

Recycling

Fragile: handle with care

Temperature limit: -20 °C ~ +55 °C

Waste electrical and electronic equipment

Serial number

Authorised representative in the European Union

Refer to the user manual/instruction booklet

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**Guilin Mecco Medical Instrument Co.,Ltd.**
**Share info GmbH**